

OPEN STAPLING. REVOLUTIONIZED.

Ten years of proven clinical performance in minimally invasive surgery, now on your open stapler

GIA™ Stapler
with Tri-Staple™ Technology



Medtronic
Further, Together

RELIABLE RESULTS. BECAUSE IT'S PROVEN TECHNOLOGY.

STRENGTH

Compared to a two-row linear stapler, the GIA™ stapler with Tri-Staple™ technology provides stronger staple line leak pressure.^{1,†}

CONFIDENCE

To optimize performance in thick-tissue applications, the GIA™ stapler with Tri-Staple™ technology offers²:

- A new knife blade with every firing
- Less retraction force than two-row linear staplers^{†‡}

FLEXIBILITY

The GIA™ stapler with Tri-Staple™ technology can be used with interchangeable cartridges to accommodate different tissue applications²

ORDERING INFORMATION

GIA™ Stapler with Tri-Staple™ Technology

	DESCRIPTION	PURPLE (MEDIUM/THICK)	BLACK (EXTRA THICK)
60 mm	Stapler	GIA60MTS ^Ω	GIA60XTS ^Ω
	Cartridge	GIA60MTC ^Ω	GIA60XTC ^Ω
80 mm	Stapler	GIA80MTS	GIA80XTS
	Cartridge	GIA80MTC	GIA80XTC

^ΩCurrently, not available for sale in the United States.



LESS STRESS

On tissue during compression and clamping^{3,‡}



GREATER PERFUSION

May be allowed into the staple line^{4,5,†,§}



CONSISTENT PERFORMANCE

Over a broad range of tissue thicknesses^{1,†}

For more information, contact your local Medtronic sales representative today.

[†]Compared to GIA™ staplers with DST Series™ technology (two-row linear staplers). [‡]Compared to GIA™ staplers with DST Series™ technology and Ethicon™ linear cutter and proximate linear cutter. [§]Preclinical results may not correlate with clinical performance in humans.

1. Based on internal test report #RE00171002 rev 0, 80 mm GIA™ stapler with Tri-Staple™ technology purple and black design verification report. Aug. 14, 2019. **2.** Based on internal test report #RE00218526, R&D memo. Aug. 27, 2019. **3.** Based on report # RE00231875, Lily tissue compression comparison. Dec. 4, 2019. **4.** Based on internal test report #RE00222215, Endo GIA™ with Tri-Staple™ technology testings performed and design similarities between Lily and Endo GIA™ with Tri-Staple™ technology memo. Oct. 16, 2019. **5.** Eschbach M, Sindberg GM, Godek ML, et al. Micro-CT imaging as a method for comparing perfusion in graduated-height and single-height surgical staple lines. *Med Devices (Auckl)*. 2018;11:267–273.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further. Together are trademarks of Medtronic. [™]* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 04/2020 – US200236 – [WF#3505420]